

Product Datasheet

Rat MOG ELISA Kit (Ready to Use) (orb555265)

Description	Rat Myelin Oligodendrocyte Glycoprotein (MOG) Sandwich ELISA Kit Ready to Use: Pre-diluted detection ...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q63345
Sensitivity	0.057ng/mL
Entrez	24558
Expiration Date	Please enquire.