

Product Datasheet

Rat NT-ProANP ELISA Kit (Ready to Use) (orb555279)

Description	Rat N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.148ng/mL
Expiration Date	Please enquire.