

Product Datasheet

Rat NPY5R ELISA Kit (Ready to Use) (orb555308)

Description	Rat Neuropeptide Y Receptor Y5 (NPY5R) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q63634
Sensitivity	0.051ng/mL
Entrez	25340
Expiration Date	Please enquire.