

Product Datasheet

Rat PTHrP ELISA Kit (Ready to Use) (orb555361)

Description	Rat Parathyroid Hormone Related Protein (PTHrP) Sandwich ELISA Kit Ready to Use: Pre-diluted detectio...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P13085
Sensitivity	4.54pg/mL
Entrez	24695
Expiration Date	Please enquire.