

Product Datasheet

Rat PGA ELISA Kit (Ready to Use) (orb555365)

Description	Rat Pepsinogen A (PGA) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient...
Reactivity	Rat
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9JJX2
Sensitivity	15.5pg/mL
Entrez	60372
Expiration Date	Please enquire.