

Product Datasheet

Rat PPAR γ C1 α ELISA Kit (Ready to Use) (orb555375)

Description	Rat Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha (PPAR γ C1 α) Sandwich ELISA K...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9QYK2
Sensitivity	0.061ng/mL
Entrez	83516
Expiration Date	Please enquire.