

Product Datasheet

Rat PLA2G2A ELISA Kit (Ready to Use) (orb555389)

Description	Rat Phospholipase A2, Group IIA (PLA2G2A) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P14423
Sensitivity	0.116ng/mL
Entrez	29692
Expiration Date	Please enquire.