

Product Datasheet

Rat PIP ELISA Kit (Ready to Use) (orb555421)

Description	Rat Prolactin Induced Protein (PIP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents f...
Reactivity	Rat
Range	46.8-3000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O70417
Sensitivity	16.23pg/mL
Entrez	64673
Expiration Date	Please enquire.