

Product Datasheet

Rat PPP2R4 ELISA Kit (Ready to Use) (orb555444)

Description	Rat Protein Phosphatase 2A Activator, Regulatory Subunit 4 (PPP2R4) Sandwich ELISA KitReady to Use: ...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	B2RYQ2
Sensitivity	0.057ng/mL
Entrez	362102
Expiration Date	Please enquire.