

Product Datasheet

Rat RLN ELISA Kit (Ready to Use) (orb555470)

Description	Rat Relaxin (RLN) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient and ...
Reactivity	Rat
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01347
Sensitivity	9.16pg/mL
Entrez	25616
Expiration Date	Please enquire.