

Product Datasheet

Rat S100A8 ELISA Kit (Ready to Use) (orb555488)

| | |
|--------------------------|---|
| Description | Rat S100 Calcium Binding Protein A8 (S100A8) Sandwich ELISA KitReady to Use: Pre-diluted detection r... |
| Reactivity | Rat |
| Range | 0.937-60ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Uniprot ID | P50115 |
| Sensitivity | 0.39ng/mL |
| Entrez | 116547 |
| Expiration Date | Please enquire. |