

Product Datasheet

Rat SCGN ELISA Kit (Ready to Use) (orb555494)

Description	Rat Secretagogen (SCGN) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien...
Reactivity	Rat
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6R556
Sensitivity	24.7pg/mL
Entrez	306942
Expiration Date	Please enquire.