

## Product Datasheet

### Rat SFRP4 ELISA Kit (Ready to Use) (orb555497)

<b>Description</b>	Rat Secreted Frizzled Related Protein 4 (SFRP4) Sandwich ELISA Kit Ready to Use: Pre-diluted detectio...
<b>Reactivity</b>	Rat
<b>Range</b>	0.187-12ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9JLS4</a>
<b>Sensitivity</b>	0.06ng/mL
<b>Entrez</b>	<a href="#">89803</a>
<b>Expiration Date</b>	Please enquire.