

Product Datasheet

Rat Serotonin Transporter ELISA Kit (Ready to Use) (orb555508)

Description	Rat Serotonin Transporter (Serotonin Transporter) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P31652
Sensitivity	0.055ng/mL
Entrez	25553
Expiration Date	Please enquire.