

Product Datasheet

Rat SDHB ELISA Kit (Ready to Use) (orb555551)

Description	Rat Succinate Dehydrogenase Complex Subunit B (SDHB) Sandwich ELISA KitReady to Use: Pre-diluted det...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.041ng/mL
Expiration Date	Please enquire.