

Product Datasheet

Rat CD71 ELISA Kit (Ready to Use) (orb555600)

Description	Rat Transferrin Receptor (CD71) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for c...
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q99376
Sensitivity	0.236ng/mL
Entrez	64678
Expiration Date	Please enquire.