

Product Datasheet

Rat TRPM4 ELISA Kit (Ready to Use) (orb555607)

Description	Rat Transient Receptor Potential Cation Channel Subfamily M, Member 4 (TRPM4) Sandwich ELISA KitRead...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9ESQ5
Sensitivity	0.055ng/mL
Entrez	171143
Expiration Date	Please enquire.