

## Product Datasheet

### Rat TFF2 ELISA Kit (Ready to Use) (orb555614)

<b>Description</b>	Rat Trefoil Factor 2 (TFF2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
<b>Reactivity</b>	Rat
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q09030</a>
<b>Sensitivity</b>	0.064ng/mL
<b>Entrez</b>	<a href="#">116592</a>
<b>Expiration Date</b>	Please enquire.