

Product Datasheet

Rat PGP9.5 ELISA Kit (Ready to Use) (orb555639)

Description	Rat Ubiquitin Carboxyl Terminal Hydrolase L1 (PGP9.5) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Rat
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q00981
Sensitivity	25pg/mL
Entrez	29545
Expiration Date	Please enquire.