

Product Datasheet

Rat UCN3 ELISA Kit (Ready to Use) (orb555647)

Description	Rat Urocortin 3 (UCN3) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	A1YKY5
Sensitivity	0.056ng/mL
Entrez	498791
Expiration Date	Please enquire.