

Product Datasheet

Rat VCC1 ELISA Kit (Ready to Use) (orb555667)

Description	Rat VEGF Co Regulated Chemokine 1 (VCC1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reage...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D4A875
Sensitivity	5.7pg/mL
Entrez	308436
Expiration Date	Please enquire.