

Product Datasheet

Primate CRTAM protein (orb555741)

Description	Recombinant Cynomolgus monkey Cytotoxic and regulatory T-cell molecule is produced by our Mammalian ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Uniprot ID	A0A2K5TKL4
Expiration Date	6 months from date of receipt.