

Product Datasheet

Primate Colony CD3E protein (orb555758)

Description	Recombinant Cynomolgus monkey T-cell surface glycoprotein CD3 epsilon chain is produced by our Mamma...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50 mM Tris-HCl, 100 mM Glycine, pH 7.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Uniprot ID	Q95LI5
Expiration Date	6 months from date of receipt.