

Product Datasheet

Mouse VEGF164 protein (orb555837)

Description	Recombinant Mouse Vascular Endothelial Growth Factor A is produced by our Yeast expression system an...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB,250mM NaCl.pH7.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Uniprot ID	Q00731
Expiration Date	6 months from date of receipt.