

Product Datasheet

GPR164 antibody (orb555911)

Description

Rabbit polyclonal antibody to GPR164

Species/Host

Rabbit

Reactivity

Human, Porcine

Conjugation

Unconjugated

Tested

IHC, IHC-P, WB

Applications

Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the center region of human GPR164. The exact sequence is proprietary.

Target

GPR164

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 50% Glycerol

Form/Appearance

Liquid: 1XPBS, 50% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Raji

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

35

Uniprot ID

[Q8TCB6](#)

Dilution Range


WB: 1:500-1:3000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of OVCA xen...



Western blot analysis of Raji cells usin...