

Product Datasheet

SEC13L1 antibody (orb555937)

Description

Rabbit polyclonal antibody to SEC13L1

Image not found or type unknown

Species/Host

Rabbit

Immunofluorescence
analysis of HepG2
cel...

Reactivity

Human

Conjugation

Unconjugated


Tested

Applications

ICC, IHC-P, IP, WB

Immunogen

Recombinant protein encompassing a sequence within the center region of human SEC13L1. The exact sequence is proprietary.



Immunohistochemical
staining of A549
xen...

Target

SEC13L1

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance


Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.



Immunoprecipitation
analysis of HepG2
wh...

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

36

Uniprot ID

P55735

Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000 , IP: 1:100-1:1000

Expiration Date

12 months from date of receipt.