

Product Datasheet

UQCRFS1 antibody (orb556021)

Description

Rabbit polyclonal antibody to UQCRFS1


Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ICC, WB

Applications

Immunofluorescence analysis of H1299 cell...

Immunogen

Recombinant protein encompassing a sequence within the center region of human UQCRFS1. The exact sequence is proprietary.

Target


UQCRFS1

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.



Immunofluorescence analysis of A549 cell...

Concentration

0.32 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Mouse brain , Raji

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


30

Uniprot ID
[P47985](https://www.uniprot.org/uniprot/P47985)
Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis of Raji cells usin...