


Product Datasheet

Rho/Rac guanine nucleotide exchange factor 18 Antibody (orb556029)

Description

Rabbit polyclonal antibody to ARHGEF18



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications



Immunofluorescence analysis of A431 cell...

Immunogen


Recombinant protein encompassing a sequence within the center region of human ARHGEF18. The exact sequence is proprietary.

Target

ARHGEF18

Preservatives

0.01% Thimerosal, 1XPBS (pH7), 1% rAlbumin, 20% Glycerol



Immunohistochemical staining of rat sple...

Form/Appearance


Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.43 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.



Immunohistochemical staining of mouse th...

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

131

Uniprot ID

Q6ZSZ5

Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.