


Product Datasheet

L3MBTL histone methyl-lysine binding protein 1 Antibody (orb556116)

Description

Rabbit polyclonal antibody to L3MBTL1



Immunofluorescence analysis of HeLa cell...

Species/Host

Rabbit

Reactivity

Human, Other

Conjugation

Unconjugated

Tested

ICC, IHC-Fr, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human L3MBTL1. The exact sequence is proprietary.

Target

L3MBTL1

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: HeLaS3, human L3MBTL-transfected 293T cells

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

84

Uniprot ID


[Q9Y468](https://www.uniprot.org/entry/Q9Y468)

Dilution Range


WB: 1:500-1:3000, IF: 1:100-1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis of 293T whole cell...



Western blot analysis of HeLa S3 cells u...