

Product Datasheet

YKL-39 antibody (orb556144)

Description

Rabbit polyclonal antibody to YKL-39

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the center region of human YKL-39. The exact sequence is proprietary.

Target

YKL-39

Preservatives

0.01% Thimerosal. 1XPBS, 1% rAlbumin, 20% Glycerol (pH7)

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: human CHI3L2-transfected 293T cells

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


44

Uniprot ID
Q15782
Dilution Range


WB: 1:1000-1:10000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of SW480
Xe...



Western blot analysis
of 293T cells usin...