

Product Datasheet

**steroidogenic acute regulatory protein Antibody
(orb556275)**

Description

Rabbit polyclonal antibody to StAR


Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications

Immunofluorescence analysis of HepG2 cel...

Immunogen

Recombinant protein encompassing a sequence within the center region of human StAR. The exact sequence is proprietary.

Target


StAR

Preservatives

0.025% ProClin 300, 1XPBS pH 7, 1% rAlbumin, 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.



Immunohistochemical staining of human he...

Concentration

0.69 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: STAR-transfected 293T

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


32

Uniprot ID
[P49675](https://www.uniprot.org/protein/P49675)
Dilution Range

WB: 1:100-1:10000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis of Mouse tissue us...