
Product Datasheet

RNA polymerase III subunit A Antibody (orb556318)

Description

Rabbit polyclonal antibody to POLR3A

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ICC, IP, WB

Immunogen

Recombinant protein encompassing a sequence within the N-terminus region of human POLR3A. The exact sequence is proprietary.

Target

POLR3A

Preservatives

0.025% ProClin 300, 1XPBS (pH7), 1% rAlbumin, 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.

Concentration

1.31 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


156

Uniprot ID
[O14802](https://www.uniprot.org/entry/O14802)
Dilution Range


WB: 1:500-1:3000, IF: 1:100-1:1000, IP: 1:100-1:500

Expiration Date

12 months from date of receipt.




Immunofluorescence analysis of HeLa cell...



Immunoprecipitation analysis of SK-N-SH

...



Western blot analysis of Various whole c...