

Product Datasheet

ZNF398 antibody (orb556417)

Description

Rabbit polyclonal antibody to ZNF398

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human ZNF398. The exact sequence is proprietary.

Target

ZNF398

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.81 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: MCF-7 , MDA-MB-231 , mouse kidney , rat kidney

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

71

Uniprot ID


Q8TD17

Dilution Range


WB: 1:500-1:3000 , IF: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of Rat tissue usin...



Western blot analysis of Mouse tissue us...