

Product Datasheet

profilin 2 Antibody (orb556545)

Description

Rabbit polyclonal antibody to Profilin 2

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human Profilin 2. The exact sequence is proprietary.

Target

Profilin 2

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: A431 , H1299 , Mouse brain , rat brain

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

15

Uniprot ID


[P35080](#)

Dilution Range


WB: 1:500-1:3000 , IHC-P: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human COX2



Western blot analysis of A431 (Lane 1), H1299 (Lane 2), Mouse brain (Lane 3), and Rat brain (Lane 4)



Western blot analysis of A431 (Lane 1), H1299 (Lane 2), Mouse brain (Lane 3), and Rat brain (Lane 4)