



Product Datasheet

phosphoserine phosphatase Antibody (orb556588)

Description	Rabbit polyclonal antibody to PSPH
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Full length human PSPH Recombinant protein.
Target	PSPH
Preservatives	0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol
Form/Appearance	Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration	1 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
MW	25
Uniprot ID	P78330
Dilution Range	WB: 1:500-1:3000 , IF: 1:100-1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of A431 cells usin...