

## Product Datasheet

### glutathione S-transferase theta 1 Antibody (orb556596)

## Description

Rabbit polyclonal antibody to GSTT1



## Species/Host

Rabbit

## Reactivity

Human, Mouse, Rat

## Conjugation

Unconjugated

## Tested

IHC-P, IP, WB

## Applications

## Immunogen

Full length human GSTT1 Recombinant protein.

## Target

GSTT1

## Preservatives

0.025% ProClin 300. 1XPBS (pH7), 20% Glycerol

## Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

## Concentration

1 mg/ml

## Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

## Note

For research use only

## Application notes

Positive control: rat liver

## Isotype

IgG

## Clonality

Polyclonal

## Purity

Purified by antigen-affinity chromatography.

## MW

27

## Uniprot ID


[P30711](#)

## Dilution Range


WB: 1:500-1:3000 , IHC-P: 1:100-1:1000 , IP: 1:100-1:500

## Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Immunoprecipitation analysis of HepG2 whole cell lysate.



Western blot analysis of mouse spleen lymphocytes.