


Product Datasheet

peroxisomal biogenesis factor 6 Antibody (orb556639)

Description	Rabbit polyclonal antibody to PEX6	
Species/Host	Rabbit	
Reactivity	Human	Western blot analysis of A549 cells using...
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	Recombinant protein encompassing a sequence within the center region of human PEX6. The exact sequence is proprietary.	
Target	PEX6	
Preservatives	0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol	
Form/Appearance	Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.	
Concentration	1 mg/ml	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.	
Note	For research use only	
Application notes	Positive control: A549	
Isotype	IgG	
Clonality	Polyclonal	
Purity	Purified by antigen-affinity chromatography.	
MW	104	
Uniprot ID	Q13608	
Dilution Range	WB: 1:500-1:3000	
Expiration Date	12 months from date of receipt.	