
Product Datasheet

SART1 antibody (orb556731)

Description

Rabbit polyclonal antibody to SART1

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the C-terminus region of human SART1. The exact sequence is proprietary.

Target

SART1

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: H1299 , HeLa , Molt-4

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


90

Uniprot ID
[O43290](#)
Dilution Range

WB: 1:1000-1:10000

Expiration Date

12 months from date of receipt.

 kDa
170
95
72
55
43


Western blot analysis of H1299(Lane 1), ...