

Product Datasheet

SSRP1 antibody (orb556766)

Description

Rabbit polyclonal antibody to SSRP1

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, IP, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human SSRP1. The exact sequence is proprietary.

Target

SSRP1

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Molt-4 , 293T , A431 , HeLa , HepG2 , Jurkat

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

81

Uniprot ID


[Q08945](https://www.uniprot.org/entry/Q08945)

Dilution Range


WB: 1:500-1:3000 , IHC-P: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of Cal27 Xe...



Western blot analysis of Whole cell extr...



Western blot analysis of Various whole c...