
Product Datasheet

sestrin 3 Antibody (orb556798)

Description

Rabbit polyclonal antibody to SESN3


Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ICC, WB

Immunogen

Recombinant protein encompassing a sequence within the center region of human SESN3. The exact sequence is proprietary.

Target

SESN3

Preservatives

0.025% ProClin 300, 1XPBS pH 7, 20% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Rat2

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


57

Uniprot ID
[P58005](https://www.uniprot.org/uniprot/P58005)
Dilution Range


WB: 1:500-1:3000, IF: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Western blot analysis of Raw264.7 whole ...



Western blot analysis of HeLa cells usin...