



Product Datasheet

Factor IX antibody (orb556879)


Description	Rabbit polyclonal antibody to Factor IX
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC-P, WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human Factor IX. The exact sequence is proprietary.
Target	Factor IX
Preservatives	0.025% ProClin 300, 1XPBS (pH7), 20% Glycerol
Form/Appearance	Liquid: 1XPBS, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.
Concentration	0.95 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
MW	52
Uniprot ID	P00740
Dilution Range	IF: 1:100-1:1000 , IHC-P: 1:100-1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...



Immunohistochemical staining of mouse pa...



Immunohistochemical staining of rat live...