


Product Datasheet

RHG09 antibody (orb55710)

Description	Rabbit polyclonal antibody to RHG09.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human RHG09
Target	ARHGAP9
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	80 kDa
Uniprot ID	Q9BRR9
NCBI	NP_115885
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
RHG09,
Sample...