

Product Datasheet

Cytomegalovirus p65 antibody (orb557158)

Description	Mouse monoclonal antibody to Cytomegalovirus p65
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1-K
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.