



## Product Datasheet

### HPV-16 E6 antibody (orb557164)

<b>Description</b>	Mouse monoclonal antibody to HPV-16 E6
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC-Fr
<b>Immunogen</b>	Please contact us for more information.
<b>Preservatives</b>	PBS with 0.02% sodium azide
<b>Concentration</b>	100 µg/ml
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1-K
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	,
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence analysis of HPV-16 E6...



Immunohistochemical staining of HPV-16 E6...