



Product Datasheet

Tuberin (phospho-Tyr1571) antibody (orb5582)

Description	Tuberin (phospho-Tyr1571) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Tuberin around the phosphorylation site of Thr1571 (R(p-Y)TE)
Target	Phospho-Tuberin (Tyr1571)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	180-220 kDa
Uniprot ID	P49815
Dilution Range	ICC=1:100, Flow-Cyt=2ug/Test, WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Flow cytometry of HeLa(black line).Prima...



ICC analysis of HeLa cell using Phospho-...