
Product Datasheet

Cynomolgus TIM-3 / HAVCR2 Protein (orb55833)

Description

Cynomolgus TIM-3, Fc Tag is expressed

contains AA Ser 22 - Arg 201...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TIM-3

PreservativesTris with Glycine, Arginine and NaCl,
pH7.5**Form/Appearance**

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 45.8 kDa. The protein migrates as 50-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

EHH54703.1

Purity

85%

Source

Cynomolgus TIM-3, Fc Tag (orb55833) is expressed from human 293 cells (HEK293). It contains AA Ser 22 - Arg 201 (Accession # EHH54703.1).

MW

45.8 kDa

Uniprot ID[G7P6Q7](#)**NCBI**[EHH54703.1](#)**Expiration Date**

6 months from date of receipt.