

## Product Datasheet

### Human VSIG8 Protein (orb55839)

<b>Description</b>	Human VSIG8, His Tag is expressed from human 293 cells (HEK293). It contains AA Val 22 - Gly 263 (Ac...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	VSIG8
<b>Preservatives</b>	20 mM NaAC, 150 mM NaCl, pH4.0
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 29.0 kDa. The protein migrates as 30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NP_001013683.1
<b>Purity</b>	95%
<b>Source</b>	Human VSIG8, His Tag (orb55839) is expressed from human 293 cells (HEK293). It contains AA Val 22 - Gly 263 (Accession # NP_001013683.1).
<b>MW</b>	29.0 kDa
<b>Uniprot ID</b>	<a href="#">Q5VU13</a>
<b>NCBI</b>	<a href="#">NP_001013683.1</a>
<b>Expiration Date</b>	6 months from date of receipt.