

Product Datasheet

NKX2-5 antibody (orb55852)

Description

NKX2-5 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human NKX2-5. Please contact us for the exact immunogen sequence. The peptide is available as orb393985.

Target

NKX2-5

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

53 kDa

Uniprot ID

[P52952](#)

NCBI

[006761](#), [006770](#)

Entrez



[8685](#)

Dilution Range


WB: 1-200-2000, IHC-P: 1:100-500

Expiration Date

12 months from date of receipt.




Immunohistochemical staining of paraffin...



IF analysis of rat stomach tissue using

...



IHC-P staining of rat heart tissue using...