



Product Datasheet

HRH4 antibody (orb55853)


Description	HRH4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human HRH4. Please contact us for the exact immunogen sequence. The peptide is available as orb393986.
Target	HRH4
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	53 kDa
Uniprot ID	Q9H3N8
NCBI	006761, 006770
Entrez	8685
Dilution Range	WB: 1-200-2000, IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



WB analysis of rat bone marrow (lane 1),...



IHC-P staining of mouse bone marrow tiss...



WB analysis of human lung cancer tissue ...