

## Product Datasheet

### PLK1 (phospho-Thr210) antibody (orb559133)

<b>Description</b>	PLK1 (phospho-Thr210) antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC
<b>Immunogen</b>	KLH conjugated synthesised phosphopeptide derived from human PLK1 around the phosphorylation site of Thr210 (KK(p-T)L)
<b>Target</b>	Phospho-PLK1 (Thr210)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Recombinant
<b>MW</b>	68 kDa
<b>CAS Number</b>	Serine/threonine-protein kinase PLK1
<b>Uniprot ID</b>	<b>P53350</b>
<b>Dilution Range</b>	Flow-Cyt=1:50, WB=1:500-1000, IHC-P=1:100-500, IHC-F=1:10-50, ICC=1:10-50, IF=1:10-50
<b>Expiration Date</b>	12 months from date of receipt.